

Product datasheet for **RC205119A1V**

Human RPL6 (NM_001024662) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RPL6 (NM_001024662) AAV Particle
Tag:	Myc-DDK
Symbol:	RPL6
Synonyms:	L6; SHUJUN-2; TAXREB107; TXREB1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205119 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGGTGAAAAAGTTGAGAAGCCAGATACTAAAGAGAAGAGACCCGAAGCCAAGAAGGTTGATGCTG
GTGGCAAGGTGAAAAAGGTAACCTCAAAGCTAAAAAGCCCAAGAAGGGGAAGCCCCATTGCAGCCGCAA
CCCTGTCCTTGTGAGAGGAGTTGGCAGGTATCCCGATCTGCCATGTATTCCAGAAAGGCCATGTACAAG
AGGAAGTACTCAGCCGCTAAATCCAAGTTGAAAAAGAAAGAGGAGAAGGTTCTCGCAACTGTTACAA
AACCAGTTGGTGGTGACAAGAACGGCGGTACCCGGGTGGTTAACTTCGCAAAATGCCTAGATATTATCC
TACTGAAGATGTGCCTCGAAAGCTGTTGAGCCACGGCAAAAAACCCCTCAGTCAGCACGTGAGAAAACTG
CGAGCCAGCATTACCCCGGGACCATTCTGATCATCCTCACTGGACGCCACAGGGGCAAGAGGGTGGTTT
TCCTGAAGCAGCTGGCTAGTGGCTTATTACTTGTGACTGGACCTCTGGTCTCAATCGAGTTCCTCTACG
AAGAACACACCAGAAATTTGTCATTGCCACTTCAACCAAAATCGATATCAGCAATGTAAAAATCCCAAAA
CATCTTACTGATGCTTACTTCAAGAAGAAGAAGCTGCGGAAGCCAGACACCAGGAAGGTGAGATCTTCG
ACACAGAAAAAGAGAAATATGAGATTACGGAGCAGCGCAAGATTGATCAGAAAGCTGTGGACTCACAAT
TTTACCAAAAATCAAAGCTATTCCTCAGCTCCAGGGCTACCTGCGATCTGTGTTTGTCTGACGAATGGA
ATTTATCCTCACAATTTGGTGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205119 protein sequence Red=Cloning site Green=Tags(s) MAGEKVEKPDTEKRPEAKKVDAGGKVKKGNLKAKKPKKGKPHCSRNPVLVRGVGRYSRSAMYSRKAMYK RKYSAASKVEKKKKEKVLATVTKPVGGDKNGGTRVVKLRKMPTYPTEDVPRKLLSHGKKPFSQHVRKL RASITPGTILIIITGRHRGKRVVFLKQLASGLLLVTGPLVLNRVPLRRTHQKFVIATSTKIDISNVKIPK HLTDAYFKKKKLRKPRHQEGEIFDTEKEKYEITEQRKIDQKAVDSQILPKIKAIPQLQGYLRSVFALTNG IYPHKLVF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001024662
ORF Size:	864 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001024662.1</u>
RefSeq Size:	1148 bp
RefSeq ORF:	867 bp
Locus ID:	6128
UniProt ID:	<u>Q02878</u>
Cytogenetics:	12q24.13
MW:	32.7 kDa