

Product datasheet for **RC205124A1V**

Human PHPT1 (NM_014172) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PHPT1 (NM_014172) AAV Particle
Tag: Myc-DDK
Symbol: PHPT1
Synonyms: CGI-202; HEL-S-132P; HSPC141; PHP; PHP14
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205124 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGTGGCGGACCTCGCTCTATTCTGATGTGGACATCGACTCCGACGGCGTCTTCAAGTATGTGC
TGATCCGAGTCCACTCGGCTCCCCGCTCCGGGCTCCGGCTGCAGAGAGCAAGGAGATCGTGCGCGGCTA
CAAGTGGGCTGAGTACCATGCGGACATCTACGACAAAGTGTGGGCGACATGCAGAAGCAAGGCTGCGAC
TGTGAGTGTCTGGGCGGGCGCATCTCCACCAGAGTCAGGACAAGAAGATTCACGTGTACGGCTATT
CCATGGCCTATGGTCTGCCAGCAGCCATTTCAACTGAGAAAATCAAAGCCAAGTACCCCGACTACGA
GGTCACCTGGGCTAACGACGGCTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205124 protein sequence
Red=Cloning site Green=Tags(s)

MAVADLALIPDVIDSDGVFKYVLRVHSAPRSGAPAAESKEIVRGYKWAHEYHADIYDKVSGDMQKQGC
CECLGGGRISHQSQDKKIHVYGYSMAYGPAQHAISTEKIKAKYPDYEVTWANDGY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_014172



[View online »](#)

ORF Size:	375 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014172.2</u>
RefSeq Size:	1199 bp
RefSeq ORF:	378 bp
Locus ID:	29085
UniProt ID:	<u>Q9NRX4</u>
Cytogenetics:	9q34.3
MW:	13.8 kDa