

Product datasheet for RC205128A1V

Human DNAJB9 (NM_012328) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DNAJB9 (NM_012328) AAV Particle
Tag: Myc-DDK
Symbol: DNAJB9
Synonyms: ERdj4; MDG-1; MDG1; MST049; MSTP049
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205128 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTACTCCCCAGTCAATTTTCATCTTTGCAATCTGCATTTTAATGATAACAGAATTAATTCTGGCCT
 CAAAAAGCTACTATGATATCTTAGGTGTGCCAAATCGGCATCAGAGCGCCAAATCAAGAAGGCCTTTCA
 CAAGTTGGCCATGAAGTACCACCTGACAAAAATAAGAGCCCGGATGCTGAAGCAAAATTCAGAGAGATT
 GCAGAAGCATATGAAACACTCTCAGATGCTAATAGACGAAAAGAGTATGATACACTTGGACACAGTGCTT
 TTAGTGTGGTAAAGGACAAAGAGGTAGTGGAAGTTCTTTGAGCAGTCATTTAACTTCAATTTTGATGA
 CTTATTTAAAGACTTTGGCTTTTGGTCAAAACCAAAACACTGGATCCAAGAAGCGTTTTGAAAATCAT
 TTCCAGACACGCCAGGATGGTGGTTCAGTAGACAAAGGCATCATTTCCAAGAATTTCTTTTGAGGTG
 GATTATTTGATGACATGTTTGAAGATATGGAGAAAATGTTTCTTTTAGTGGTTTGACTCTACCAATCA
 GCATACAGTACAGACTGAAAATAGATTTTCATGGATCTAGCAAGCACTGCAGGACTGTCACTCAACGAAGA
 GGAAATATGGTTACTACATACACTGACTGTTCAGGACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205128 protein sequence
 Red=Cloning site Green=Tags(s)

MATPQSIFIFAICILMITELILASKSYDILGVPKSASERQIKKAFHKLAMKYHPDNKNSPDAAKFREI
 AEAYETLSDANRRKEYDTLGHSAFTSGKGQRGSGSSFEQSFNFDDLFKDFGFFGQNQNTGSKKRFENH
 FQTRQDGGSSRRHHFQEFSGGGFLDDMFEDMEKMFSGFDSTNQHTVQTENRFHGSSKHCRTVTQRR
 GNMVTTYTDCSQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_012328
ORF Size:	669 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_012328.1
RefSeq Size:	2538 bp
RefSeq ORF:	672 bp
Locus ID:	4189
UniProt ID:	Q9UBS3
Cytogenetics:	14q24.2-q24.3
MW:	25.5 kDa