

## Product datasheet for **RC205130L3V**

### **Cbl c (CBLC) (NM\_012116) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cbl c (CBLC) (NM_012116) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cbl c
Synonyms:	CBL-3; CBL-SL; RNF57
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012116
ORF Size:	1422 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205130).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012116.2</a>
RefSeq Size:	1591 bp
RefSeq ORF:	1425 bp
Locus ID:	23624
UniProt ID:	<a href="#">Q9ULV8</a>
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome



[View online »](#)

<b>Protein Pathways:</b>	Chronic myeloid leukemia, Endocytosis, ErbB signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, Pathways in cancer, T cell receptor signaling pathway, Ubiquitin mediated proteolysis
<b>MW:</b>	52.5 kDa
<b>Gene Summary:</b>	This gene encodes a member of the Cbl family of E3 ubiquitin ligases. Cbl proteins play important roles in cell signaling through the ubiquitination and subsequent downregulation of tyrosine kinases. Expression of this gene may be restricted to epithelial cells, and alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2012]