

Product datasheet for **RC205148L4V**

Brevican (BCAN) (NM_021948) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Brevican (BCAN) (NM_021948) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Brevican |
| Synonyms: | BEHAB; CSPG7 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_021948 |
| ORF Size: | 2733 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC205148). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_021948.3 |
| RefSeq Size: | 3470 bp |
| RefSeq ORF: | 2736 bp |
| Locus ID: | 63827 |
| UniProt ID: | Q96GW7 |
| Cytogenetics: | 1q23.1 |
| Domains: | CCP, Xlink, CLECT, ig, IGv, IG, EGF, EGF |
| Protein Families: | Secreted Protein |


[View online »](#)

MW: 99.1 kDa

Gene Summary: This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]