

## Product datasheet for **RC205159L3V**

### **TBC1D16 (NM\_019020) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TBC1D16 (NM_019020) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TBC1D16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019020
ORF Size:	2301 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205159).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019020.2</a>
RefSeq Size:	10956 bp
RefSeq ORF:	2304 bp
Locus ID:	125058
UniProt ID:	<a href="#">Q8TBP0</a>
Cytogenetics:	17q25.3
MW:	86.4 kDa
Gene Summary:	May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]


[View online »](#)