

Product datasheet for **RC205171L4V**

TMEM90B (SYNDIG1) (NM_024893) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM90B (SYNDIG1) (NM_024893) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM90B
Synonyms:	C20orf39; DSPC2; IFITMD5; TMEM90B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024893
ORF Size:	774 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024893.1
RefSeq Size:	2523 bp
RefSeq ORF:	777 bp
Locus ID:	79953
UniProt ID:	Q9H7V2
Cytogenetics:	20p11.21
Protein Families:	Transmembrane
MW:	28.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that belongs to the interferon-induced transmembrane family of proteins. A similar protein in rat is thought to regulate the development of excitatory synapses. [provided by RefSeq, Jul 2013]