

Product datasheet for **RC205205L3V**

PHYHIP (NM_014759) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PHYHIP (NM_014759) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHYHIP
Synonyms:	DYRK1AP3; PAHX-AP; PAHXAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014759
ORF Size:	990 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205205).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014759.2 , NP_055574.2
RefSeq Size:	3232 bp
RefSeq ORF:	993 bp
Locus ID:	9796
UniProt ID:	Q92561
Cytogenetics:	8p21.3
Domains:	FN3
MW:	37.6 kDa


[View online »](#)

Gene Summary:

Its interaction with PHYH suggests a role in the development of the central system.
[UniProtKB/Swiss-Prot Function]