

Product datasheet for **RC205208L4V**

KIFAP3 (NM_014970) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KIFAP3 (NM_014970) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIFAP3
Synonyms:	dj190I16.1; FLA3; KAP-1; KAP-3; KAP3; SMAP; Smg-GDS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014970
ORF Size:	2376 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205208).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014970.3
RefSeq Size:	3011 bp
RefSeq ORF:	2379 bp
Locus ID:	22920
UniProt ID:	Q92845
Cytogenetics:	1q24.2
Domains:	Armadillo_seg
Protein Families:	Druggable Genome



[View online »](#)

MW: 91.2 kDa

Gene Summary: The small G protein GDP dissociation stimulator (smg GDS) is a regulator protein having two activities on a group of small G proteins including the Rho and Rap1 family members and Ki-Ras; one is to stimulate their GDP/GTP exchange reactions, and the other is to inhibit their interactions with membranes. The protein encoded by this gene contains 9 'Armadillo' repeats and interacts with the smg GDS protein through these repeats. This protein, which is highly concentrated around the endoplasmic reticulum, is phosphorylated by v-src, and this phosphorylation reduces the affinity of the protein for smg GDS. It is thought that this protein serves as a linker between human chromosome-associated polypeptide (HCAP) and KIF3A/B, a kinesin superfamily protein in the nucleus, and that it plays a role in the interaction of chromosomes with an ATPase motor protein. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]