

Product datasheet for **RC205211L4V**

THEG (NM_016585) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	THEG (NM_016585) Human Tagged ORF Clone Lentiviral Particle
Symbol:	THEG
Synonyms:	CT56; THEG1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016585
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205211).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016585.3
RefSeq Size:	1341 bp
RefSeq ORF:	1140 bp
Locus ID:	51298
UniProt ID:	Q9P2T0
Cytogenetics:	19p13.3
MW:	43.4 kDa


[View online »](#)

Gene Summary:

This gene is specifically expressed in the nucleus of haploid male germ cells. The orthologous gene in mice encodes a protein that may play a role in protein assembly through interactions with T-complex protein 1 subunit epsilon. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]