

Product datasheet for RC205215A1V

Human C20orf141 (NM_080739) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C20orf141 (NM_080739) AAV Particle
Tag: Myc-DDK
Symbol: C20orf141
Synonyms: dj860F19.4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205215 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACCCGGCTCTGCTTACCCAGACCCGAAGCACGTGAGGATCCGATCCAGTTCCTCCAAGGGGCCTGG
 GTGCTGGGGAGGGGTCAAGTAGTCCAGTGCGTCCACCTGTATCCACCTGGGGCCCTAGCTGGGCCAGCT
 CCTGGACAGTGTCTATGGCTGGGGGCACTAGGACTGACAATCCAGGCAGTCTTTCCACCACTGGCCCA
 GCCCTGCTGCTGCTTCTGGTCACTTCTCACCTTTGACCTGCTCCATAGGCCCGCAGGTCACTCTGC
 CACAGCGCAAACCTTCTACCAGGGGCCAGAGTCAGGGGGCCGGTGAAGTCTGGACAGCAGGAGGCTCT
 ACTCCTGCAAATGGGTACAGTCTCAGGACAACCTAGCCTCCAGGACGCACTGCTGCTGCTCATGGGG
 CTGGGCCCGCTCTGAGAGCCTGTGGCATGCCCTTGACCTGCTTGGCCTGGCTTTCTGCCTCCATCCTT
 GGGCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205215 protein sequence
 Red=Cloning site Green=Tags(s)

MTRLCLPRPEAREDPVPPRGLGAGEGSGSPVPPVSTWGPSWAQLLDVSLWLGALGLTIQAVFSTTGP
 ALLLLLVSFLTFDLLHRPAGHTLPQRKLLTRGQSQGAGEGPGQGEALLLQMGTVSGQLSLQDALLLLLMG
 LGPLLACGMPLTLLGLAFCLHPWA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_080739
ORF Size:	495 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080739.2</u> , <u>NP_542777.1</u>
RefSeq Size:	741 bp
RefSeq ORF:	498 bp
Locus ID:	128653
UniProt ID:	<u>Q9NUB4</u>
Cytogenetics:	20p13
MW:	17.4 kDa