

## Product datasheet for **RC205225L2V**

### OGDHL (NM\_018245) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	OGDHL (NM_018245) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OGDHL
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018245
ORF Size:	3030 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205225).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018245.1</a>
RefSeq Size:	3774 bp
RefSeq ORF:	3033 bp
Locus ID:	55753
UniProt ID:	<a href="#">Q9ULD0</a>
Cytogenetics:	10q11.23
Protein Pathways:	Citrate cycle (TCA cycle), Lysine degradation, Metabolic pathways, Tryptophan metabolism
MW:	114.5 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is similar to oxoglutarate dehydrogenase (OGDH) of the OGDH complex, which degrades glucose and glutamate. This gene encodes several isoforms, including some that appear to localize to mitochondria. The encoded protein down-regulates the AKT signaling cascade and can suppress the growth of cervical cancer cells. [provided by RefSeq, Dec 2016]