

Product datasheet for **RC205228L4V**

KLHL11 (NM_018143) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLHL11 (NM_018143) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLHL11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018143
ORF Size:	2124 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018143.1
RefSeq Size:	2326 bp
RefSeq ORF:	2127 bp
Locus ID:	55175
UniProt ID:	Q9NVR0
Cytogenetics:	17q21.2
Domains:	BTB, Kelch
Protein Families:	Secreted Protein
MW:	80.1 kDa



[View online »](#)

Gene Summary:

Component of a cullin-RING-based BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex that mediates the ubiquitination of target proteins, leading most often to their proteasomal degradation.[UniProtKB/Swiss-Prot Function]