

Product datasheet for **RC205229L4V**

ZNF334 (NM_199441) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF334 (NM_199441) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF334
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199441
ORF Size:	1926 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205229).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199441.1 , NP_955473.1
RefSeq Size:	4932 bp
RefSeq ORF:	1929 bp
Locus ID:	55713
Cytogenetics:	20q13.12
Protein Families:	Transcription Factors
MW:	75 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the C2H2 zinc finger family. The encoded protein contains a Krueppel-associated box, fourteen C2H2 zinc finger domains, and four C2H2-type/integrase DNA-binding domains. Decreased expression of this gene may be a marker for rheumatoid arthritis. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jul 2012]