

Product datasheet for **RC205251L4V**

DNAJC5B (NM_033105) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC5B (NM_033105) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC5B
Synonyms:	CSP-beta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033105
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205251).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033105.4
RefSeq Size:	1398 bp
RefSeq ORF:	600 bp
Locus ID:	85479
UniProt ID:	Q9UF47
Cytogenetics:	8q13.1
MW:	22.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the DNAJ heat shock protein 40 family of co-chaperone proteins that is characterized by an N-terminal DNAJ domain, a linker region, and a cysteine-rich C-terminal domain. The encoded protein, together with heat shock protein 70, is thought to regulate the proper folding of other proteins. The orthologous mouse protein is membrane-associated and is targeted to the trans-golgi network. [provided by RefSeq, Mar 2017]