

Product datasheet for **RC205258A1V**

Human CRHBP (NM_001882) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CRHBP (NM_001882) AAV Particle
Tag:	Myc-DDK
Symbol:	CRHBP
Synonyms:	CRF-BP; CRFBP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205258 representing NM_001882 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCGCCCAACTTCAAACCTTCAGTGTCACTTCATTCTCATCTTCCTGACGGCTCTAAGAGGGGAAAGCC
GGTACCTAGAGCTGAGGGAAGCGGCGGACTACGATCCTTTCTGCTCTTCAGCGCCAACCTGAAGCGGGA
GCTGGCTGGGAGCAGCCGTACCGCCGCGCTCTGCGGTGCCTGGACATGCTGAGCCTCCAGGGCCAGTTC
ACCTTCACCGCCGACCGGCCGAGCTGCACTGCGCAGCCTTCTTCATCAGCGAGCCCGAGGAGTTCATTA
CCATCCACTACGACCAAGGTCTCCATCGACTGTCAGGGCGGCGACTTCTGAAGGTATTTGATGGTTGGAT
TCTCAAGGGGGAGAAGTTCCCCAGTTCCAGGATCATCTCTCCCTCAGCTGAGCGGTACATAGATTTTC
TGTGAGAGTGGTCTTAGCAGGAGGAGCATCAGATCTTCCAGAATGTGGCCATGATCTTCTTCCGAGTCC
ATGAACCAGGAAATGGATTACATTAACCATAAAGACAGACCCCAACCTTTCTTGAATGTCATTTTC
TCAGACTCCAAATGGAAAGTTTACCCTGGTAGTTCCACACCAGCATCGAACTGCAGCTTCTCCATAATT
TATCCTGTGGTGATCAAATATCTGATCTTACCCTGGGACACGTAATGGTCTTCAGTTAAAGAAATCCT
CAGCAGGTTGCGAGGGAATAGGAGACTTTGTGGAGCTGCTGGGAGGAAGTGGATTGGACCCTTCCAAGAT
GACGCCTTTAGCTGATCTCTGCTACCCCTTTCATGGCCCGGCCAGATGAAAGTTGGCTGTGACAACACT
GTGGTGCGCATGGTCTCCAGTGGAAAACACGTAATCGTGTGACTTTTGAGTATCGTCAGCTGGAGCCGT
ACGAGCTGGAACCCAAATGGAACAGTATCGGGGAATTCTGTTTGTCTGGTCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205258 representing NM_001882 Red=Cloning site Green=Tags(s) MSPNFKLQCHFILIFLTALRGESRYLELREAADYDPFLLFSANLKRELAGEQPYRRALRCLDMLSLQGQF TFTADRPQLHCAAFFISEPEEFITIHVDQVSIDCQGGDFLKVFDGWILKGEKFPSSQDHPLPSAERYIDF CESGLSRRSIRSSQNVAMIFFRVHEPGNGFTLTIKTDPNLFPCNVISQTPNGKFTLVVPHQHRNCSFSII YPVVIKISDLTLGHVNGQLKKSSAGCEGIGDFVELLGGTGLDPSKMTPLADLCYPFHGPAQMKVGCNT VVRMVSSGKHVNRVTFEYRQLEPYELENPNNGNSIGEFCLSGL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001882
ORF Size:	966 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001882.3</u>
RefSeq Size:	1854 bp
RefSeq ORF:	969 bp
Locus ID:	1393
UniProt ID:	<u>P24387</u>
Cytogenetics:	5q13.3
MW:	36 kDa