

Product datasheet for **RC205263L1V**

Protein phosphatase inhibitor 1 (PPP1R1A) (NM_006741) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Protein phosphatase inhibitor 1 (PPP1R1A) (NM_006741) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Protein phosphatase inhibitor 1 |
| Synonyms: | I1; IPP1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_006741 |
| ORF Size: | 513 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC205263). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006741.2 , NP_006732.2 |
| RefSeq Size: | 1849 bp |
| RefSeq ORF: | 516 bp |
| Locus ID: | 5502 |
| UniProt ID: | Q13522 |
| Cytogenetics: | 12q13.2 |
| Protein Families: | Druggable Genome |



[View online »](#)

Protein Pathways: Long-term potentiation

MW: 18.9 kDa

Gene Summary: Inhibitor of protein-phosphatase 1. This protein may be important in hormonal control of glycogen metabolism. Hormones that elevate intracellular cAMP increase I-1 activity in many tissues. I-1 activation may impose cAMP control over proteins that are not directly phosphorylated by PKA. Following a rise in intracellular calcium, I-1 is inactivated by calcineurin (or PP2B). Does not inhibit type-2 phosphatases.[UniProtKB/Swiss-Prot Function]