

Product datasheet for **RC205271A1V**

Human RASL10A (NM_006477) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RASL10A (NM_006477) AAV Particle
Tag: Myc-DDK
Symbol: RASL10A
Synonyms: RRP22
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205271 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGGGTAGCCTGCGGGTGGCCGTTCTAGGCGCCCCGGGCGTGGGCAAGACGGCCATCATCCGCCAGT
TCGTGTTTCGGTGACTACCCCGAGCGCCACCGGCCACGGACGGGCCGCGCCTCTACCGACCCGCGGTGCT
GCTCGACGGCGCCGCTCTACGACTTGAGCATCCGCGACGGCGACGTCGCTGGCCCCGGCTCGAGCCCCGGG
GGTCCGGAGGAGTGGCCAGACGCTAAGGACTGGAGCTTGACAGGACACGGACGCCCTTCGTGCTCGTCTACG
ACATCTGCAGCCCGACAGTTTCGACTACGTGAAGGCCCTGCGGCAGCGCATCGCGGAGACCAGGCCGGC
GGGCGCGCCCGAAGCGCCCATCCTCGTGGTAGGCAACAAGCGGGACAGGACAGCGGCTGCGCTTCGGACCG
CGGCGCGCGCTGGCCGCCCTAGTGCGCAGGGGCTGGCGCTGCGGCTACCTCGAGTGCTCCGCCAAGTACA
ACTGGCACGTGCTGCGTCTCTCCGCGAGCTGCTGCGCTGCGCTCTGGTGCAGCGCGCCCTGCACACCC
GGCCCTGCGCTGCAGGGGGCGCTGCATCCCGCGCGCTGCAGCCTCATG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205271 protein sequence
Red=Cloning site Green=Tags(s)

MGGLSLRVAVLGAPGVGKTAIRQFVFGDYPERHRPTDGPRLYRPAVLLDGAVYDLSIRDGDVAGPGSSPG
GPEEWPDAKDWSLQDTDAFVLVDICSPDSFDYVKALRQRIAE TRPAGAPEAPILVVGNGRDRQRLRFGP
RRALAAVRRGWRCGYLECSAKYNWHVLRFLRELLRCALVRARPAHPALRLQALHPARCSLM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_006477
ORF Size:	609 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006477.3
RefSeq Size:	1504 bp
RefSeq ORF:	612 bp
Locus ID:	10633
UniProt ID:	Q92737
Cytogenetics:	22q12.2
MW:	22.5 kDa