

Product datasheet for **RC205271L3V**

RASL10A (NM_006477) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RASL10A (NM_006477) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RASL10A
Synonyms:	RRP22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006477
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205271).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006477.3
RefSeq Size:	1504 bp
RefSeq ORF:	612 bp
Locus ID:	10633
UniProt ID:	Q92737
Cytogenetics:	22q12.2
Domains:	RAS, RHO, RAB
MW:	22.5 kDa



[View online »](#)

Gene Summary: Potent inhibitor of cellular proliferation.[UniProtKB/Swiss-Prot Function]