

## Product datasheet for **RC205282L3V**

### EDIL3 (NM\_005711) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EDIL3 (NM_005711) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EDIL3
Synonyms:	DEL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005711
ORF Size:	1440 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205282).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005711.3</a>
RefSeq Size:	4771 bp
RefSeq ORF:	1443 bp
Locus ID:	10085
UniProt ID:	<a href="#">O43854</a>
Cytogenetics:	5q14.3
Domains:	F5_F8_type_C, EGF_CA, EGF, EGF
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**MW:** 53.8 kDa

**Gene Summary:** The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior. [provided by RefSeq, Jul 2008]