

Product datasheet for **RC205293A1V**

Human NEUROD2 (NM_006160) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NEUROD2 (NM_006160) AAV Particle
Tag:	Myc-DDK
Symbol:	NEUROD2
Synonyms:	bHLHa1; DEE72; EIEE72; NDRF
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205293 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTGACCCGCTGTTTCAGCGAGCCCGCCTTCTCTCGGACGTGCCCAAGTTCGCCAGCTGGGGCGACG
GCGAAGACGACGAGCCGAGGAGCGACAAGGGCGACGCGCCGCCACCGCCACCGCTGCGCCCGGGCCAGG
GGCTCCGGGGCCAGCCCGGGCGGCCAAGCCAGTCCCTCTCCGTGGAGAAGAGGGGACGGAGGCCACGTTG
GCCGAGGTCAAGGAGGAAGGCGAGCTGGGGGGAGAGGAGGAGGAGGAAGAGGAGGAGGAAGAAGGACTGG
ACGAGGCGGAGGGCGAGCGGCCAAGAAGGGCGGCCAAGAAGCGCAAGATGACCAAGCGCGCTTGGA
GCGCTCCAAGCTTCGGCGGCAGAAAGGCGAACGCGGGAGCGCAACCGCATGCACGACCTGAACGCAGCC
CTGGACAACCTGCGCAAGGTGGTGCCCTGCTACTCCAAGACGCGAGAAGCTGTCCAAGATCGAGACGCTGC
GCCTAGCCAAGAATACTATCTGGGCGCTCTCGGAGATCCTGCGCTCCGGCAAGCGGCCAGACCTAGTGTC
CTACGTGCAGACTCTGTGCAAGGTCTGTGCGAGCCACCAACCTGCTGGTGGCCGGCTGTCTGCAGCTC
AACTCTCGCAACTTCTCACGGAGCAAGGCGCCGACGGTGCCGGCCGCTTCCACGGCTCGGGCGGCCGT
TCGCCATGCACCCCTACCGTACCGTGCTGCGCCTGGCGGGCGCACAGTGCCAGGCGGCCGGCGCCCT
GGGCGGCGGCGGCGCACGCCCTGCGACCCACGGCTACTGCGCAGCCTACGAGACGCTGTATGCGCG
GCAGGCGGTGGCGGCGGAGCCCGACTACAACAGCTCCGAGTACGAGGGCCCGCTCAGCCCCCGCTCT
GTCTCAATGGCAACTTCTCACTCAAGCAGGACTCCTCGCCGACCACGAGAAAAGCTACCACTACTCTAT
GCACTACTCGGCGCTGCCGGTTTCGCGGCCACGGGCCACGGGCTAGTCTTCGGCTCGTGGCTGTGCGC
GGGGCGTCCACTCGGAGAATCTTTGTCTTACGATATGCACCTTCACCACGACCGGGGGCCCATGTACG
AGGAGCTCAATGCGTTTTTTCATAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205293 protein sequence Red=Cloning site Green=Tags(s) MLTRLFSEPGLLSDVPKFASWGDGEDDEPRSDKGDAPPPPPAPGPGAPGPARAAKPVPLRGEETATL AEVKEEGELGGEEEEEEEEGLDEAGERPKKGGPKRKMTKARLERSKLRRQKANARERNRMHDLNAA LDNLRKVVPYCYSKTQKLSKIETLRLAKNYIWALSEILRSGKRPDLVSYVQTLCKGLSQPTTNLVAGCLQL NSRNFLTEQGADGAGRFHSGGPFAMHPYPYPCSRLAGAQCQAAGGLGGGAHALRTHGYCAAYETLYAA AGGGGASPDYNSSEYEGPLSPPLCLNGNFSLKQDSSPDHEKSYHYSMHYSALPGSRPTGHGLVFGSSAVR GGVHSENLLSYDMHLHHDGPMYEELNAFFHN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006160
ORF Size:	1146 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006160.3, NP_006151.2</u>
RefSeq Size:	3048 bp
RefSeq ORF:	1149 bp
Locus ID:	4761
UniProt ID:	<u>Q15784</u>
Cytogenetics:	17q12
MW:	41.3 kDa