

## Product datasheet for **RC205314L4V**

### **RIOK3 (NM\_145906) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RIOK3 (NM_145906) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RIOK3
Synonyms:	SUDD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145906
ORF Size:	1272 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205314).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145906.1</a> , <a href="#">NP_665913.1</a>
RefSeq Size:	3589 bp
RefSeq ORF:	1274 bp
Locus ID:	8780
Cytogenetics:	18q11.2
Domains:	RIO
Protein Families:	Druggable Genome, Protein Kinase
MW:	48.4 kDa



[View online »](#)

**Gene Summary:**

This gene was first identified by the similarity of its product to the *Aspergillus nidulans* Sudd protein. This gene is now recognized as a member of the right open reading frame (RIO) kinase gene family. This gene encodes a serine/threonine kinase that localizes to the cytoplasm and plays a role in the processing of the pre-40 S ribosomal subunit. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2017]