

## Product datasheet for **RC205324L3V**

### ID2 (NM\_002166) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ID2 (NM_002166) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ID2
Synonyms:	bHLHb26; GIG8; ID2A; ID2H
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002166
ORF Size:	402 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205324).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002166.4</a>
RefSeq Size:	1402 bp
RefSeq ORF:	405 bp
Locus ID:	3398
UniProt ID:	<a href="#">Q02363</a>
Cytogenetics:	2p25.1
Domains:	HLH



[View online »](#)

<b>Protein Families:</b>	ES Cell Differentiation/IPS, Stem cell relevant signaling - TGFb/BMP signaling pathway, Transcription Factors
<b>Protein Pathways:</b>	TGF-beta signaling pathway
<b>MW:</b>	14.9 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the inhibitor of DNA binding family, members of which are transcriptional regulators that contain a helix-loop-helix (HLH) domain but not a basic domain. Members of the inhibitor of DNA binding family inhibit the functions of basic helix-loop-helix transcription factors in a dominant-negative manner by suppressing their heterodimerization partners through the HLH domains. This protein may play a role in negatively regulating cell differentiation. A pseudogene of this gene is located on chromosome 3. [provided by RefSeq, Aug 2011]