

Product datasheet for **RC205372L3V**

PYK2 (PTK2B) (NM_173174) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PYK2 (PTK2B) (NM_173174) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PYK2 |
| Synonyms: | CADTK; CAKB; FADK2; FAK2; PKB; PTK; PYK2; RAFTK |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_173174 |
| ORF Size: | 3027 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC205372). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_173174.1 |
| RefSeq Size: | 4715 bp |
| RefSeq ORF: | 3030 bp |
| Locus ID: | 2185 |
| UniProt ID: | Q14289 |
| Cytogenetics: | 8p21.2 |
| Protein Families: | Druggable Genome, Protein Kinase |



[View online »](#)

| | |
|--------------------------|---|
| Protein Pathways: | Calcium signaling pathway, Chemokine signaling pathway, GnRH signaling pathway, Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity |
| MW: | 115.7 kDa |
| Gene Summary: | <p>This gene encodes a cytoplasmic protein tyrosine kinase which is involved in calcium-induced regulation of ion channels and activation of the map kinase signaling pathway. The encoded protein may represent an important signaling intermediate between neuropeptide-activated receptors or neurotransmitters that increase calcium flux and the downstream signals that regulate neuronal activity. The encoded protein undergoes rapid tyrosine phosphorylation and activation in response to increases in the intracellular calcium concentration, nicotinic acetylcholine receptor activation, membrane depolarization, or protein kinase C activation. This protein has been shown to bind CRK-associated substrate, nephrocystin, GTPase regulator associated with FAK, and the SH2 domain of GRB2. The encoded protein is a member of the FAK subfamily of protein tyrosine kinases but lacks significant sequence similarity to kinases from other subfamilies. Four transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> |