

## Product datasheet for **RC205375L4V**

### AGFG1 (NM\_004504) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | AGFG1 (NM_004504) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | AGFG1  |
| Synonyms:                 | HRB; RAB; RIP  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_004504  |
| ORF Size:                 | 1686 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC205375).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_004504.3</a>  |
| RefSeq Size:              | 8684 bp  |
| RefSeq ORF:               | 1689 bp  |
| Locus ID:                 | 3267   |
| UniProt ID:               | <a href="#">P52594</a>   |
| Cytogenetics:             | 2q36.3   |
| Domains:                  | ArfGap   |
| MW:                       | 58.3 kDa   |



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008]