

## Product datasheet for **RC205391L2V**

### **ADH1B (NM\_000668) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ADH1B (NM_000668) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADH1B
Synonyms:	ADH2; HEL-S-117
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000668
ORF Size:	1125 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205391).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000668.3</a>
RefSeq Size:	2707 bp
RefSeq ORF:	1128 bp
Locus ID:	125
UniProt ID:	<a href="#">P00325</a>
Cytogenetics:	4q23
Domains:	ADH_zinc_N
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Drug metabolism - cytochrome P450, Fatty acid metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Retinol metabolism, Tyrosine metabolism

**MW:** 39.9 kDa

**Gene Summary:** The protein encoded by this gene is a member of the alcohol dehydrogenase family. Members of this enzyme family metabolize a wide variety of substrates, including ethanol, retinol, other aliphatic alcohols, hydroxysteroids, and lipid peroxidation products. This encoded protein, consisting of several homo- and heterodimers of alpha, beta, and gamma subunits, exhibits high activity for ethanol oxidation and plays a major role in ethanol catabolism. Three genes encoding alpha, beta and gamma subunits are tandemly organized in a genomic segment as a gene cluster. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2013]