

Product datasheet for **RC205422A1V**

Human TEX47 (NM_152706) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TEX47 (NM_152706) AAV Particle
Tag:	Myc-DDK
Symbol:	TEX47
Synonyms:	C7orf62
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205422 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCCTTTTCGGTCCATAACCAGAAGGGCAGCAAAAGGCCTTTGCCACTGGAACCTCTTCTTTTCTCC
AAGTCCCACGTAGCAATTACCTGCACCTTTCAAGAAGAGAAACAACGACTACACCTAAAGAAATTCCTTCT
TGATAGGATGTTTCTAGTGGCCAAGATACAAGCAATGTAGAAAGAAAAGATGTTGCTGACTACTATGAA
CAAATGTTTCAGTCAGTTTTGAAACATCACCTAGGAGAAGCAGTGACAGGATTGCTGCTCATCTATCCCA
CTTCCATTCTGCATATCCTCGAGTCCTCCAGCGACACTCTCTACAAAGTTCTTTTAGATTATATTGGCCA
TGTCAAAGATGAAACAGTATTTTTTATTCAACAAATGAAAATTATAGTCATTTCTATAACATTCCAATG
AGGCTTTTTATGCAATGGCATGTTTCAGTGATAAAAGTTCCAGTTATGTATCTCGACGATGTGACACAGT
CACAGTCCCTAAAGGAGGTATCACAGATTTTCTCACACAACTCATAAACTGTCACTCTACCTTTGCCA
GACTATGAAAGTAGGCACTAAAGGACCAGGCGATAACTTACACCAAGTTGCACCTGACCTACTCTCCCA
GAACAAATCATAAAGTACTTGTGCAATCCGAAGAATTCATGGACCCGGCAACATTTATAAATGTATA
ATAGACCCATACATCACTCTGGATTCTGAGGTGGTATGGCCTGCTCCTTCACGTTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC205422 protein sequence Red=Cloning site Green=Tags(s) MSFSVHNQKGSKRPLPLELLFLQVPRSNYLHFQEEKQRLHLKKFLLDRMFLVAKIQANVERKDVADYYE QMFQSVLKHHLGEAVTGLLLIYPTSLIHILESSDTLYKVLDDYIGHVKDETVFFIQQMKIIVISHNIPM RLFMQWHVSVIKVPVMYLDDVTQSQLKEVITDFLTQTHKLSLYLCQTMKVGTKGPGDNLHQVAPDLLLP EQIIKYLCCKSEEFMDPATFINMYNRPIHITLDSEVVWPAPSRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152706
ORF Size:	759 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152706.2</u>
RefSeq Size:	1038 bp
RefSeq ORF:	762 bp
Locus ID:	219557
UniProt ID:	<u>Q8TBZ9</u>
Cytogenetics:	7q21.13
MW:	29.5 kDa