

Product datasheet for **RC205424L3V**

SYNPR (NM_144642) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SYNPR (NM_144642) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SYNPR
Synonyms:	SPO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144642
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205424).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144642.3
RefSeq Size:	2632 bp
RefSeq ORF:	798 bp
Locus ID:	132204
UniProt ID:	Q8TBG9
Cytogenetics:	3p14.2
Domains:	MARVEL
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

MW: 29.2 kDa

Gene Summary: Intrinsic membrane protein of small synaptic vesicles. Probable vesicular channel protein (By similarity).[UniProtKB/Swiss-Prot Function]