

Product datasheet for **RC205428A1V**

Human KHDRBS2 (NM_152688) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human KHDRBS2 (NM_152688) AAV Particle
Tag:	Myc-DDK
Symbol:	KHDRBS2
Synonyms:	SLM-1; SLM1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205428 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGAAGAGGAGAAATATTTGCCTGAGCTGATGGCAGAGAAAGATAGCCTGGATCCATCTTTTGTGCATG
CGTCGCGCCTTTTGGCAGAAGAAATTGAAAAGTTTCAAGGTTCTGATGGAAAAAGGAAGACGAAGAAAA
GAAGTATCTTGATGTCATCAGCAACAAAAACATAAAGCTCTCAGAAAGAGTACTGATTCCTGTCAAGCAG
TATCCAAAGTTCAATTTTGTGGGAAATTGCTTGGACCAAGAGGAACTCCTTGAAGAGGCTACAGGAAG
AAACAGGTGCTAAAATGTCTATCCTGGGCAAAGGATCAATGAGAGATAAAGCTAAGGAAGAAGAACTAAG
GAAGAGTGGGGAAGCCAAATATGCCCACTTGAGTGATGAGCTTCATGTATTAATTGAAGTGTGTTGCTCCA
CCTGGGGAAGCTTATTCACGTATGAGTCATGCATTGGAAGAGATTAAAAAATTCCTGGTTCCTGACTACA
ATGATGAAATTCGTGAGGAACAACACTACGTGAATTATCTTAAATGGCTCAGAGGACTCTGGTCGTGG
CAGAGGTATTAGAGGCAGAGGGATCAGAATAGCTCCACAGCTCCTTCAAGGGGCCGTGGGGTGCCATT
CCTCTCCCCCACCACCTGGACGAGGTGTTCTACCCCTCGGGGAAGCACTGTAACCCGTGGAGCGCTTC
CAGTGCCACCTGTAGCAAGAGGTGTCCCTACCCCTCGAGCCCGGGGGGCACCAACAGTGCCAGGATACAG
GGCACCTCCTCCTCCAGCCCATGAAGCTTATGAAGAATATGGTTATGATGATGGCTACGGGGGTGAATAT
GATGACCAGACCTATGAGACTTATGATAACAGCTACGCGACCCAAACACAAAGTGTGCCTGAATACTATG
ACTACGGTCATGGAGTAAGTGAGGATGCCTATGACAGCTACGCACCAGAAGAATGGGCCACAACCAGCTC
TAGCTTGAAGGCACCACCGCAAAGGTACGCCAGAGGGGATACAGGGAACACCCCTATGTTAGATAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205428 protein sequence Red=Cloning site Green=Tags(s) MEEEKYLPELMAEKDSLDPFVHASRLLAEEIEKFQGS DGKKEDEEKKYLDVISNKNIKLSERVLIPVKQ YPKFNFVGKLLGPRGNSLKRLQEETGAKMSILGKGS MRDKAKEEELRKS GEAKY A H L S D E L H V L I E V F A P PGEAYSRMSHALEEIKKFLVPDYNDEIRQEQLRELSYLN GSEDSGRGRGIRGRGIRIAPTAPSRGRGGAI PPPPPPGRGVLTPRGSTVTRGALPVPPVARGVPTPRARGAPTVPGYRAPPPPAHEAYEEYGYDDGYGGEY DDQTYETYDNSYATQTQSVPEYYDYGHGVSEDAYDSYAPEEWATTSSSLKAPPQRSARGGYREHPYGRY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152688
ORF Size:	1047 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152688.1, NP_689901.1</u>
RefSeq Size:	2336 bp
RefSeq ORF:	1050 bp
Locus ID:	202559
UniProt ID:	<u>Q5VWX1</u>
Cytogenetics:	6q11.1
MW:	38.9 kDa