

## Product datasheet for **RC205435L4V**

### **DPP10 (NM\_020868) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DPP10 (NM_020868) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DPP10
Synonyms:	DPL2; DPPY; DPRP-3; DPRP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020868
ORF Size:	2388 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205435).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020868.2</a>
RefSeq Size:	5253 bp
RefSeq ORF:	2391 bp
Locus ID:	57628
UniProt ID:	<a href="#">Q8N608</a>
Cytogenetics:	2q14.1
Domains:	Peptidase_S9, DPPIV_N_term
Protein Families:	Druggable Genome, Protease, Transmembrane


[View online »](#)

**MW:** 90.9 kDa

**Gene Summary:** This gene encodes a single-pass type II membrane protein that is a member of the S9B family in clan SC of the serine proteases. This protein has no detectable protease activity, most likely due to the absence of the conserved serine residue normally present in the catalytic domain of serine proteases. However, it does bind specific voltage-gated potassium channels and alters their expression and biophysical properties. Mutations in this gene have been associated with asthma. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]