

## Product datasheet for **RC205438L4V**

### Prion protein PrP (PRNP) (NM\_000311) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Prion protein PrP (PRNP) (NM_000311) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Prion protein PrP
Synonyms:	AltPrP; ASCR; CD230; CJD; GSS; KURU; p27-30; PRIP; PrP; PrP27-30; PrP33-35C; PrPc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000311
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205438).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000311.3</a>
RefSeq Size:	2755 bp
RefSeq ORF:	762 bp
Locus ID:	5621
UniProt ID:	<a href="#">P04156</a>
Cytogenetics:	20p13
Domains:	prion
Protein Families:	ES Cell Differentiation/IPS, Stem cell - Pluripotency, Transmembrane



[View online »](#)

**Protein Pathways:** Prion diseases

**MW:** 27.6 kDa

**Gene Summary:** The protein encoded by this gene is a membrane glycosylphosphatidylinositol-anchored glycoprotein that tends to aggregate into rod-like structures. The encoded protein contains a highly unstable region of five tandem octapeptide repeats. This gene is found on chromosome 20, approximately 20 kbp upstream of a gene which encodes a biochemically and structurally similar protein to the one encoded by this gene. Mutations in the repeat region as well as elsewhere in this gene have been associated with Creutzfeldt-Jakob disease, fatal familial insomnia, Gerstmann-Straussler disease, Huntington disease-like 1, and kuru. An overlapping open reading frame has been found for this gene that encodes a smaller, structurally unrelated protein, AltPrp. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2014]