

Product datasheet for **RC205471L4V**

Tumor protein D52 like 3 (TPD52L3) (NM_001001875) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tumor protein D52 like 3 (TPD52L3) (NM_001001875) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Tumor protein D52 like 3
Synonyms:	D55; hD55; NYDSP25; TPD55
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001001875
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205471).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001001875.2 , NP_001001875.1
RefSeq Size:	1353 bp
RefSeq ORF:	399 bp
Locus ID:	89882
UniProt ID:	Q96J77
Cytogenetics:	9p24.1
Protein Families:	Druggable Genome



[View online »](#)

MW: 14.6 kDa

Gene Summary: This gene encodes a member of the tumor protein D52-like family of proteins. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. The encoded protein may play a role in spermatogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]