

Product datasheet for **RC205472A1V**

Human SRRP35 (SRSF12) (NM_080743) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SRRP35 (SRSF12) (NM_080743) AAV Particle
Tag:	Myc-DDK
Symbol:	SRRP35
Synonyms:	SFRS13B; SFRS19; SRrp35
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205472 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCCTCGCTACACGAGGCCCAACACCTCCCTGTTTCATCAGGAACGTCGCGGACGCCACCAGGCCTG
AGGACTTGCGCCGTGAGTTTGGTCGATATGGCCCTATAGTAGACGTTTACATTCCACTTGACTTCTACAC
TCGCCGCCCAAGAGGATTTGCTTATGTTCAATTTGAAGATGTTTCGAGGTGCTGAAGATGCTCTTTATAAC
CTCAATAGAAAGTGGGTATGTGGCCGTCAGATTGAAATACAGTTTGCACAAGGTGATCGAAAAACACAG
GCCAAATGAAATCAAAAGAAGTCATCCTTGTCTCCAAGTGATCACAGGAGATCAAGAAGCCCCAGCCA
AAGAAGAAGTCAAGTAGAAGTTCTTCATGGGGAAGAAATAGGAGGCGGTGACAGAGCCTTAAAGAGTCT
CGACACAGGCGATTTTCTTATAGCAAGTCTAAATCTCGTTCCAAATCATTACCAAGGCGGTCTACCTCAG
CAAGGCAGTCAAGAAGTCCAAGAAGGAATTTGGCTCTAGAGGACGGTCAAGGTCCAAGTCCCTTACAAAA
GAGGTCCAAGTCAATAGGAAATCACAGTCAAGTTCACCTCAAAGCAGACTAGCTCAGGAACAAATCA
AGATCACATGGAAGACATTCTGACTCAATAGCAAGATCCCGTGTAATCTCCAAAGGGTATACCAATT
CTGAACTAAAGTACAAACAGCAAGCATTCTCATTTTCGGTCACATTCCAGATCTCGAAGTTATCGTCA
TAAAAACAGTTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205472 protein sequence Red=Cloning site Green=Tags(s) MSRYTRPPNTSLFIRNVADATRPEDLRREFGRYGPIVDVYIPLDFYTRRPRGFAYVQFEDVRGAEDALYN LNRKWVCGRQIEIQFAQGDRKTPGQMKSKERHPCSPSDHRRSRSPSQRRTRSRSSSWGRRRRSDSLKES RHRRFYSYKSKSRKSLPRRSTSAQSRTPRRNFGSRGRSRKSLQKRKSIGKSQSSSPQKQTSSGTKS RSHGRHSDSIARSPCKSPKGYTNSETKVQTAKHSHFRSHSRSRSYRHKNSW TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_080743
ORF Size:	783 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080743.2</u>
RefSeq Size:	3628 bp
RefSeq ORF:	786 bp
Locus ID:	135295
UniProt ID:	<u>Q8WXF0</u>
Cytogenetics:	6q15
MW:	30.5 kDa