

Product datasheet for RC205494L2

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Non Neuronal Enolase (ENO1) (NM_001428) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Non Neuronal Enolase (ENO1) (NM_001428) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Non Neuronal Enolase

Synonyms: ENO1L1; HEL-S-17; MPB1; NNE; PPH

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC205494).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001428

ORF Size: 1302 bp



ORÏGENE

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001428.2</u>

 RefSeq Size:
 2204 bp

 RefSeq ORF:
 1305 bp

 Locus ID:
 2023

 UniProt ID:
 P06733

 Cytogenetics:
 1p36.23

Domains: enolase

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Glycolysis / Gluconeogenesis, Metabolic pathways, RNA degradation

MW: 47.2 kDa


Gene Summary: This gene encodes alpha-enolase, one of three enolase isoenzymes found in mammals. Each

isoenzyme is a homodimer composed of 2 alpha, 2 gamma, or 2 beta subunits, and functions as a glycolytic enzyme. Alpha-enolase in addition, functions as a structural lens protein (taucrystallin) in the monomeric form. Alternative splicing of this gene results in a shorter isoform that has been shown to bind to the c-myc promoter and function as a tumor suppressor. Several pseudogenes have been identified, including one on the long arm of chromosome 1. Alpha-enolase has also been identified as an autoantigen in Hashimoto encephalopathy.

[provided by RefSeq, Jan 2011]



Product images:



Circular map for RC205494L2