

Product datasheet for **RC205503A1V**

Human NANP (NM_152667) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NANP (NM_152667) AAV Particle
Tag:	Myc-DDK
Symbol:	NANP
Synonyms:	C20orf147; dj694B14.3; HDHD4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205503 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGCTGAGCCGCGTGCGGGCGGTTTTCTTTGACTTGGACAACACTCTCATCGACACGGCCGGGGCGA
GCAGGAGAGGCATGTTGGAGGTGATAAACTCTTACAATCAAAATACCATTATAAAGAAGAGGCTGAAAT
CATCTGTGATAAAGTTCAAGTTAACTCAGCAAGGAATGTTTTCATCCTTACAATACATGCATTACTGAT
TTAAGGACTTCACATTGGGAAGAAGCAATCCAGGAAACAAAAGGTGGTGCAGCCAATAGAAAATTGGCTG
AAGAATGTTATTTCTTTGAAATCTACAGTTTACAGCATATGACACTAGCAGAAGACGTCAAAGCCAT
GCTTACTGAACCTTCGAAAGGAGGTCCGCCTACTTCTATTAACGAATGGGGACAGACAGCCAGAGGGAG
AAGATTGAGGCTTGTGCCTGTCAGTCTATTTTACGCTGTTGTTGTAGGTGGAGAGCAGAGAGAGGAGA
AACCAGCACCGTCCATATTTTATTACTGCTGCAATCTTCTCGGAGTACAACCTGGGGACTGTGTGATGGT
CGGTGACACATTAGAAACCGACATCCAAGGAGGCCTCAATGCAGGATTGAAAGCAACAGTCTGGATCAAT
AAAAATGGAATAGTGCCACTGAAGTCTCCCGATTCCGCATTACATGGTTTCTTCTGTGCTAGAGTTAC
CTGCTCTCTTACAAAGTATAGACTGCAAAGTCAGTATGTCCACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205503 protein sequence Red=Cloning site Green=Tags(s) MGLSRVRAVFFDLNLTIDTAGASRRGMLEVIKLLQSKYHYKEEAEIICDKVQVKLSKECFHPYNTCITD LRTSHWEEAIQETKGGGAANRKLAEECYFLWKSTRLQHMTLAEDVKAMLTELKREVRLLLLTNGDRQTQRE KIEACACQSYFDAVVVGGEQREEKPAPSIFFYYCCNLLGVQPGDCVMVGDLETDIQGGLNAGLKATVWIN KNGIVPLKSSPVPHYMVSSVLELPALLQSIDCKVSMST TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152667
ORF Size:	744 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152667.1</u>
RefSeq Size:	3806 bp
RefSeq ORF:	747 bp
Locus ID:	140838
UniProt ID:	<u>Q8TBE9</u>
Cytogenetics:	20p11.21
MW:	27.8 kDa