

## Product datasheet for **RC205504L3V**

### STK11IP (NM\_052902) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	STK11IP (NM_052902) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STK11IP
Synonyms:	LIP1; LKB1IP; STK11IP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_052902
ORF Size:	3297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205504).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_052902.2</a> , <a href="#">NP_443134.1</a>
RefSeq Size:	3602 bp
RefSeq ORF:	3267 bp
Locus ID:	114790
UniProt ID:	<a href="#">Q8N1F8</a>
Cytogenetics:	2q35
Domains:	LRR
Protein Families:	Druggable Genome


[View online »](#)

**MW:** 121.3 kDa

**Gene Summary:** May regulate STK11/LKB1 function by controlling its subcellular localization.[UniProtKB/Swiss-Prot Function]