

Product datasheet for **RC205512L4V**

DIS3L2 (NM_152383) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DIS3L2 (NM_152383) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DIS3L2
Synonyms:	FAM6A; hDIS3L2; PRLMNS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152383
ORF Size:	2655 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205512).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152383.4
RefSeq Size:	3507 bp
RefSeq ORF:	2658 bp
Locus ID:	129563
UniProt ID:	Q8IYB7
Cytogenetics:	2q37.1
MW:	99.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is similar in sequence to 3'/5' exonucleolytic subunits of the RNA exosome. The exosome is a large multimeric ribonucleotide complex responsible for degrading various RNA substrates. Several transcript variants, some protein-coding and some not, have been found for this gene. [provided by RefSeq, Mar 2012]