

Product datasheet for **RC205519A1V**

Human C1orf19 (TSEN15) (NM_052965) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human C1orf19 (TSEN15) (NM_052965) AAV Particle
 Tag: Myc-DDK
 Symbol: C1orf19
 Synonyms: C1orf19; PCH2F; sen15
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC205519 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGAGCGCGCGATTCCGAGCCGACCCCGGCTGCAGCGGCCTGGGTCCGGGCGGTGTTCCGCGCT
 TTGGCGACGGCGGTGGAGCTCCTTCGTGGGCCCTGAGGACGCCTGGATGGGCACTCACCTAAGTATCT
 AGAAATGATGGAATTAGATATAGGAGATGCCACCAAGTTTATGTAGCGTTCTTGTTTACCTGGACCTC
 ATGGAAGCAAAAGCTGGCATGAAGTAACTGTGTAGGATTACCAGAACTCCAGCTCATCTGCCTTGTG
 TACTGAGATAGAAGGGGAGGGTTACAGACTGTGGTGCCTACCCCATCACTGCTTCCCTCAGCCATAA
 CAGGATAAGGGAGATCTTGAAGCATCTCGAAAGTTGCAAGGTATCCAGATTTGCCGATGTCTTTTACT
 TTGGCCATAGTGGAGTCTGATTCTACAATAGTCTATTATAAACTTACTGATGGATTTATGCTGCCAGACC
 CTCAGAATATTCTCTTAGAAGA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205519 protein sequence
 Red=Cloning site Green=Tags(s)

MEERGDSEPTPGCSGLPGGVRFVGDGGGAPSWAPEDAWMGTHPKYLEMMELDIGDATQVYVAFVYLDL
 MESKSWHEVNCVGLPELQLICLVGTEIEGELQTVVPTPITASLSHNRIREILKASRKLQGDPLPMSFT
 LAIVESDSTIVYYKLTGDFMLPDPQNISLRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
 Serotype: AAV-2



[View online »](#)

ACCN:	NM_052965
ORF Size:	513 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_052965.1</u>
RefSeq Size:	1998 bp
RefSeq ORF:	516 bp
Locus ID:	116461
UniProt ID:	<u>Q8WW01</u>
Cytogenetics:	1q25.3
MW:	18.6 kDa