

Product datasheet for **RC205527L4V**

SHCBP1 (NM_024745) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SHCBP1 (NM_024745) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SHCBP1
Synonyms:	PAL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024745
ORF Size:	2016 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205527).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024745.2
RefSeq Size:	3249 bp
RefSeq ORF:	2019 bp
Locus ID:	79801
UniProt ID:	Q8NEM2
Cytogenetics:	16q11.2
Domains:	PbH1
MW:	75.6 kDa



[View online »](#)

Gene Summary:

May play a role in signaling pathways governing cellular proliferation, cell growth and differentiation. May be a component of a novel signaling pathway downstream of Shc. Acts as a positive regulator of FGF signaling in neural progenitor cells.[UniProtKB/Swiss-Prot Function]