

Product datasheet for **RC205533A1V**

Human TMEM163 (NM_030923) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TMEM163 (NM_030923) AAV Particle
Tag:	Myc-DDK
Symbol:	TMEM163
Synonyms:	DC29; SV31
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205533 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCCGGCCGCGGGCATCCAGCGCCGCAGCTCCCAGGGGCCACCGTCCCGCCGCCGCCCGGGGCC
ACGCGCCACCGGCTGCCGCCCCGCGCCCGGCCCGCTGAGCTCCCGGTGCGCGAGCCGCCAGCTGGA
GGAGGAGCGGCAGGTGCGGATCAGCGAGAGCGGCCAGTTCAGCGACGGGCTGGAGGACCGAGGCTTACTA
GAAAGCAGCACCCGCTGAAACCTCACGAAGCCAGAACTACAGGAAGAAGGCATTGTGGGTGTCCTGGT
TCTCCATCATTGTACCCTGGCCCTCGCGGTGGCTGCCCTTACTGTCTCCGTTATGAGGTACAGCGCCTC
TGCTTTTGGGTTTGATTTGATGCCATCCTGGACGTCTGTGTCGCGGATTGTCCTGTGGCGTTACAGC
AACGCGCCGCTGTGCACTCTGCCCATAGGGAGTACATAGCCTGTGTATCTTGGGGGTGATATTCCTTC
TGTCATCCATATGTATAGTGGTCAAAGCCATCCATGACCTCTCAACTAGGCTGCTCCAGAAAGTGACGA
TTTCTGTTCAGTGTCTCCATTTAAGTGGATTCTTTGCAGCATCCTGGCCGTGTTGAAGTTCATGCTG
GGGAAGGTTCTGACCAGTAGAGCACTCATAACAGATGGGTTAACTCCCTCGTGGGTGGCGTGATGGGCT
TCTCATTCTTCTGAGCGCGAAGTGTCAAGCATGACTCAGCGGTCTGGTACCTGGACGGCAGCATAGG
CGTTCTGATCGGCCTACCATATTTGCCTATGGGTCAAACCTCCTCATCGACATGGTGCCGAGGGTGAGG
CAGACAGTCACTACGAGATGTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205533 protein sequence Red=Cloning site Green=Tags(s) MEPAAGIQRRSSQGPTVPPPPRGHAPPAAGPAPLSSPVREPPQLEERQVRISESGQFSDGLEDRGLL ESSTRLPHEAQNYRKKALWVSWFSIIIVTLALAVAAFTVSMRYASAFGFAFDAILDVLSAIVLWRY NAAVHSAHREYIACVILGVIFLLSSICIVVKAIHDLSTRLLPEVDDFLFSVILSGILCSILAVLKFML GKVLTSRALITDGFNSLVGGVMGFSILLSAEVFKHDSAVWYLDGSIGVLIGLTIFAYGVKLLIDMVPVR QTRHYEMFE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_030923
ORF Size:	867 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_030923.2</u>
RefSeq Size:	1901 bp
RefSeq ORF:	870 bp
Locus ID:	81615
UniProt ID:	<u>Q8TC26</u>
Cytogenetics:	2q21.3
MW:	31.5 kDa