

## Product datasheet for **RC205533L4V**

### TMEM163 (NM\_030923) Human Tagged ORF Clone Lentiviral Particle

#### Product data:


Product Type:	Lentiviral Particles
Product Name:	TMEM163 (NM_030923) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM163
Synonyms:	DC29; SV31
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030923
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205533).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030923.2</a>
RefSeq Size:	1901 bp
RefSeq ORF:	870 bp
Locus ID:	81615
UniProt ID:	<a href="#">Q8TC26</a>
Cytogenetics:	2q21.3
Protein Families:	Transmembrane
MW:	31.5 kDa



[View online »](#)

**Gene Summary:** May bind zinc and other divalent cations and recruit them to vesicular organelles.  
[UniProtKB/Swiss-Prot Function]

**Product images:**



[RC205533L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC205533L4V particle to overexpress human TMEM163-mGFP fusion protein.