

Product datasheet for **RC205539L4V**

GPR62 (NM_080865) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR62 (NM_080865) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR62
Synonyms:	GPCR8; KPG_005
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080865
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205539).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080865.3
RefSeq Size:	2206 bp
RefSeq ORF:	1107 bp
Locus ID:	118442
UniProt ID:	Q9BZJ7
Cytogenetics:	3p21.2
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	37.6 kDa



[View online »](#)

Gene Summary:

Orphan G-protein coupled receptor. Constitutively activates the G(q/11)/inositol phosphate and the G(s)-alpha/cAMP signaling pathways (PubMed:28827538). Has spontaneous activity for beta-arrestin recruitment (PubMed:28827538). Shows a reciprocal modulation of signaling functions with the melatonin receptor MTNR1B most likely through receptor heteromerization (PubMed:28827538).[UniProtKB/Swiss-Prot Function]