

## Product datasheet for **RC205559L2V**

### **KBTBD7 (NM\_032138) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KBTBD7 (NM_032138) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KBTBD7
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_032138
ORF Size:	2052 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205559).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032138.3</a>
RefSeq Size:	3008 bp
RefSeq ORF:	2055 bp
Locus ID:	84078
UniProt ID:	<a href="#">Q8WVZ9</a>
Cytogenetics:	13q14.11
Domains:	BTB, Kelch
MW:	77.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a transcriptional activator, having been shown to increase the transcription of activator protein-1 and serum response element. The encoded protein can also form a complex with KBTBD6 and CUL3, which regulates the ubiquitylation and degradation of TIAM1, which is a regulator of RAC1. [provided by RefSeq, Jul 2016]