

Product datasheet for **RC205591L3V**

Cyclin M3 (CNNM3) (NM_017623) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyclin M3 (CNNM3) (NM_017623) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cyclin M3
Synonyms:	ACDP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017623
ORF Size:	2121 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205591).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017623.4
RefSeq Size:	4916 bp
RefSeq ORF:	2124 bp
Locus ID:	26505
UniProt ID:	Q8NE01
Cytogenetics:	2q11.2
Protein Families:	Transcription Factors, Transmembrane
MW:	75.9 kDa



[View online »](#)

Gene Summary: Probable metal transporter.[UniProtKB/Swiss-Prot Function]