

Product datasheet for **RC205595A1V**

Human Prokineticin 1 (PROK1) (NM_032414) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human Prokineticin 1 (PROK1) (NM_032414) AAV Particle

Tag: Myc-DDK

Symbol: Prokineticin 1

Synonyms: EGVEGF; PK1; PRK1

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC205595 representing NM_032414
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGAGGTGCCACGCGAGTCTCAATCATGCTCCTCTAGTAAGTGTGTCTGACTGTGCTGTGATCACAG
GGGCCTGTGAGCGGGATGTCCAGTGTGGGCAGGCACCTGCTGTGCCATCAGCCTGTGGCTTCGAGGGCT
GCGGATGTGCACCCCGCTGGGGCGGGAAGGCGAGGAGTGCCACCCCGGCAGCCACAAGATCCCCTTCTTC
AGGAAACGCAAGCACACCTGTCTTGTCTTGCCTAACCTGCTGTGCTCCAGGTTCCCGGACGGCAGGT
ACCGTGTCTCCATGGACTTGAAGAATCAATTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205595 representing NM_032414
Red=Cloning site Green=Tags(s)

MRGATRVSIMLLLVTVSDCAVITGACERDVQCAGTCCAISLWLRGLRMCTPLGREGECHPGSHKIPFF
RKRKHHTCPCLPNLLCSRFPDGRYRCSMDLKNINF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_032414

ORF Size: 315 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032414.2</u>
RefSeq Size:	1372 bp
RefSeq ORF:	318 bp
Locus ID:	84432
UniProt ID:	<u>P58294</u>
Cytogenetics:	1p13.3
MW:	11.5 kDa