

## Product datasheet for **RC205599L2V**

### ANGEL1 (NM\_015305) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ANGEL1 (NM_015305) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANGEL1
Synonyms:	Ccr4e; KIAA0759
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015305
ORF Size:	2010 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205599).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015305.2</a> , <a href="#">NP_056120.1</a>
RefSeq Size:	4131 bp
RefSeq ORF:	2013 bp
Locus ID:	23357
UniProt ID:	<a href="#">Q9UNK9</a>
Cytogenetics:	14q24.3
Domains:	Exo_endo_phos
MW:	75.2 kDa



[View online »](#)