

Product datasheet for **RC205612L3V**

SERPINB4 (NM_002974) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SERPINB4 (NM_002974) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SERPINB4
Synonyms:	LEUPIN; PI11; SCCA-2; SCCA1; SCCA2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002974
ORF Size:	1170 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205612).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002974.2
RefSeq Size:	1787 bp
RefSeq ORF:	1173 bp
Locus ID:	6318
UniProt ID:	P48594
Cytogenetics:	18q21.33
Protein Families:	Druggable Genome, Secreted Protein
MW:	44.9 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the serpin family of serine protease inhibitors. The encoded protein is highly expressed in many tumor cells and can inactivate granzyme M, an enzyme that kills tumor cells. This protein, along with serpin B3, can be processed into smaller fragments that aggregate to form an autoantigen in psoriasis, probably by causing chronic inflammation. [provided by RefSeq, Jan 2017]