

Product datasheet for **RC205634L4V**

MS4A8B (MS4A8) (NM_031457) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | MS4A8B (MS4A8) (NM_031457) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MS4A8B |
| Synonyms: | 4SPAN4; CD20L5; MS4A4; MS4A8B |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_031457 |
| ORF Size: | 750 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC205634). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_031457.1 , NP_113645.1 |
| RefSeq Size: | 1369 bp |
| RefSeq ORF: | 753 bp |
| Locus ID: | 83661 |
| UniProt ID: | Q9BY19 |
| Cytogenetics: | 11q12.2 |
| Protein Families: | Druggable Genome, Transmembrane |
| MW: | 26.3 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the membrane-spanning 4A gene family. Members of this protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. The gene encoding this protein is localized to 11q12.3, among a cluster of family members. [provided by RefSeq, Jul 2008]