

Product datasheet for **RC205640**

ABCG2 (NM_004827) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABCG2 (NM_004827) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ABCG2
Synonyms:	ABC15; ABCP; BCRP; BCRP1; BMDP; CD338; CDw338; EST157481; GOUT1; MRX; MXR; MXR-1; MXR1; UAQTL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC205640 representing NM_004827
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCCTCCAGTAATGTGCAAGTTTTATCCCAGTGTACAAGGAAACACCAATGGCTTCCCCGCGACAG
 CTTCCAATGACCTGAAGGCATTTACTGAAGGAGCTGTGTTAAGTTTTATAACATCTGCTATCGAGTAAA
 ACTGAAGAGTGGCTTTCTACCTGTGCAAAACCAAGTTGAGAAAGAAATATTATCGAATATCAATGGGATC
 ATGAAACCTGGTCTCAACGCCATCTGGGACCCACAGGTGGAGGCAAAATCTTCGTTATTAGATGTCTTAG
 CTGCAAGGAAAGATCCAAGTGGATTATCTGGAGATGTTCTGATAAATGGAGCACCGCGACCTGCCAATTT
 CAAATGTAATTCAGGTTACGTGGTACAAGATGATGTTGTGATGGGCACTCTGACGGTGAGAGAAAACCTTA
 CAGTTCTCAGCAGCTCTTCGGCTTGCAACAATATGACGAATCATGAAAAAACGAACGGATTAAACAGGG
 TCATTCAAGAGTTAGGTCTGGATAAAGTGGCAGACTCCAAGGTTGGAACCTCAGTTTATCCGTGGTGTGTC
 TGGAGGAGAAAGAAAAGGACTAGTATAGGAATGGAGCTTACTGATCCTTCCATCTTGTCTTGGAT
 GAGCCTACAACCTGGCTTAGACTCAAGCACAGCAAAATGCTGTCCTTTGCTCCTGAAAAGGATGTCTAAGC
 AGGGACGAACAATCATCTTCTCCATTATCAGCCTCGATATTCATCTTCAAGTTGTTTGATAGCCTCAC
 CTTATTGGCCTCAGGAAGACTTATGTTCCACGGGCCTGCTCAGGAGGCCTTGGGATACTTTGAATCAGCT
 GGTTATCACTGTGAGGCCTATAATAACCTGCAGACTTCTTCTGGACATCATTAAATGGAGATTCCACTG
 CTGTGGCATTAAACAGAGAAGAAGACTTTAAAGCCACAGAGATCATAGAGCCTTCCAAGCAGGATAAGCC
 ACTCATAGAAAAATTAGCGGAGATTTATGTCAACTCCTCCTTCTACAAAGAGACAAAAGCTGAATTACAT
 CAATTTCCGGGGGTGAGAAGAAGAAGATCACAGTCTTCAAGGAGATCAGCTACACCACCTCCTTCT
 GTCATCAACTCAGATGGGTTTCCAAGCGTTCATTCAAAAACCTTGCTGGGTAATCCCCAGGCCTCTATAGC
 TCAGATCATTGTACAGTCGTACTGGGACTGGTTATAGGTGCCATTTACTTTGGGCTAAAAAATGATTCT
 ACTGGAATCCAGAACAGAGCTGGGGTTCTCTTCTCTGACGACCAACCAGTGTTTCAGCAGTGTTTCAG
 CCGTGGAATCCTTGTGGTAGAGAAGAAGCTCTTCATACATGAATACATCAGCGGATACTACAGAGTGTC
 ATCTTATTTCTTGAAAACTGTTATCTGATTTATTACCATGAGGATGTTACCAAGTATTATTTACC
 TGTATAGTGTACTTCATGTTAGGATTGAAGCCAAAGGCAGATGCCTTCTTCGTTATGATGTTTACCCTTA
 TGATGGTGGCTTATTCAGCCAGTTCATGGCACTGGCCATAGCAGCAGGTGAGAGTGTTTCTGTAGC
 AACACTTCTCATGACCATCTGTTTGTGTTTATGATGATTTTTTCAGGTCTGTTGGTCAATCTCACAAAC
 ATTGCATCTTGGCTGTCATGGCTTCAGTACTTCAGCATTCCACGATATGGATTTACGGCTTTCAGCATA
 ATGAATTTTTGGGACAAAACCTTGCCCAAGGACTCAATGCAACAGGAAACAATCCTTGTAACATATGCAAC
 ATGTACTGGCGAAGAATATTTGGTAAAGCAGGGCATCGATCTCTCACCCTGGGGCTTGTGGAAGAATCAC
 GTGGCCTTGGCTTGATGATTGTTATTTTCTCACAATTGCCTACCTGAAATTGTTATTTCTTAAAAAAT
 ATTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205640 representing NM_004827
 Red=Cloning site Green=Tags(s)


MSSSNVEVFIPVSQGNTNGFPATASNDLKAFTEGAVLSFHNICYRVKLKSGFLPCRKPVEKEILSNINGI
 MKPGLNAILGPTGGGKSSLLDVLAAKDPSSGLSGDVLINGAPRPANFKCNSGYVVQDDVVMGTLTVRENL
 QFSAALRLATTMTNHEKNERINRVIQELGLDKVADSKVGTQFIRGVSGGERKRTSIGMELITDPSILFLD
 EPTTGLDSSTANAVLLLLKRMSKQGRITIFSIIHQPRYSIFKLFDSLTLASGRLMFHGPAQEALGYFESA
 GYHCEAYNNPADFFLDIINGDSTAVLNREEDFKATEIIIEPSKQDKPLIEKLAIEIYVNSSFYKETKAELH
 QLSGGEKKKKITVFEISYTTSFCHQLRWVSKRSFKNLLGNPQASIAQIIIVTVVLGLVIGAIYFGLKNDS
 TGIQNRAGVLFLLTTNQCFSSVSAVELFVVEKKLFIHEYISGYRVSSYFLGKLLSDLLPMRMLPSIIFT
 CIVYFMLGLKPKADAFFVMMFTLMMVAYSASSMALAIAAGQSVSVATLLMTICFVMMIFSGLLVNLTT
 IASWLSWLQYFSIPRYGFTALQHNEFLGQNFPCPLNATGNNPCNYATCTGEEYLVKQGIDLSPWGLWKNH
 VALACMIVIFLTIAYLKLLFLKKYS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg2774_e01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_004827

ORF Size: 1965 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_004827.3](#)

RefSeq Size: 4445 bp

RefSeq ORF: 1968 bp

Locus ID: 9429

UniProt ID: [Q9UNQ0](#)

Cytogenetics: 4q22.1

Domains: ABC_tran, AAA


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways: ABC transporters

MW: 72.1 kDa

Gene Summary: The membrane-associated protein encoded by this gene is included in the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. Alternatively referred to as a breast cancer resistance protein, this protein functions as a xenobiotic transporter which may play a major role in multi-drug resistance. It likely serves as a cellular defense mechanism in response to mitoxantrone and anthracycline exposure. Significant expression of this protein has been observed in the placenta, which may suggest a potential role for this molecule in placenta tissue. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RC205640