

## Product datasheet for **RC205643L4V**

### CHERP (NM\_006387) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CHERP (NM_006387) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CHERP
Synonyms:	DAN16; SCAF6; SRA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006387
ORF Size:	2652 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205643).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006387.5</a> , <a href="#">NP_006378.2</a>
RefSeq Size:	4105 bp
RefSeq ORF:	2751 bp
Locus ID:	10523
UniProt ID:	<a href="#">Q8IWX8</a>
Cytogenetics:	19p13.11
Domains:	Surp, G-patch
Protein Pathways:	Spliceosome



[View online »](#)

**MW:** 100 kDa

**Gene Summary:** Involved in calcium homeostasis, growth and proliferation.[UniProtKB/Swiss-Prot Function]